

Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs

and Drug Control

MUSCAT



سِلاطِنَا جُمَاةِ
وَزَارَةِ الصِّحَّةِ
وَالرِّقَابَةِ الدَّوَلِيَّةِ
مَسْقَط

To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES

Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)

Director General of Engineering Affairs, MOH

Director General of Royal Hospital

Director General of Khoula Hospital

Director General of Medical Supplies (MOH)

Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)

Hospital Director (Al Nahda Hospital)

Hospital Director (Al Massara Hospital)

The Head of Medical Services in SQU Hospital

The Head of Medical Services in Royal Oman Police

The Head of Medical Services in Ministry of Defence

The Head of Medical Services in The Diwan

The Head of Medical Services in The Sultan's Special Force

The Head of Medical Services in Internal Security Services

The Head of Medical Services in Petroleum Development of Oman

The Head of Medical Services in LNG Oman

ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Please find attached our Circular No. 52..... dated 02/03/20 regarding recall of In-vitro diagnostics - microbiological products from Siemens Healthcare

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information

Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs
and Drug Control



سلطنة عمان
وزارة الصحة
المديرية العامة للصيد
والرقابة الدوائية
مسقط

MUSCAT

Circular No. 52 / 2020

07-07-1441 H

02-03-2020

Ref: 40/2020

Recall of In-vitro diagnostics - microbiological products from Siemens Healthcare

Source of Recall	NCMDR National Centre for Medical Devices Reporting, SFDA https://ncmdr.sfda.gov.sa/Secure/CA/CAViewRecall.aspx?caid=6&rid=15061
Product	IMMULITE 2000/IMMULITE 2000 XPi EBV-VCA IgM-In-vitro diagnostics - microbiological products
Manufacturer	Siemens Healthcare Diagnostics Product
Local Agent	New Source
The affected products	Test Code: ECM Catalog: L2KEM2 SMN: 10488005 Lot: 355
Reason for Recall	Siemens Healthcare Diagnostics Inc. has confirmed the potential for an increase in within run and within laboratory imprecision on Siemens' EBV-VCA IgM Control Module quality control material and patient samples with EBV-VCA IgM on the IMMULITE® 2000/IMMULITE® 2000 XPi Systems These samples may exhibit higher percent coefficient of variation (%CV) than the precision performance data published in the Instructions For Use across nonreactive, indeterminate and reactive S/CO ratios. The imprecision is especially evident on quality control material and samples that are close to the indeterminate range (0.9 to 1.0 S/CO ratio) of the assay. The imprecision may cause the interpretation of a sample that is nonreactive to be indeterminate or reactive and a sample that is reactive to be nonreactive or indeterminate.
Action	1. Kindly check your stock, contact your local agent for remedial action
Comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical Device to Director of Medical Device Control contact E-mail dg-padc@moh.gov.om

HN

Directorate General of Pharmaceutical Affairs & Drug Control
Sultanate of Oman



Dr. Mohammed Hamdan Al Rubaie

DIRECTOR GENERAL