

# Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs  
and Drug Control

MUSCAT



سِلاطِنَة عُمان  
وَزارة الصِّحة  
والمَدِيْرَة العَامَة للصِّدْقَة  
والمَرْقَبَة الدَّوْلِيَّة  
مَسَقَط

To:

**THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES**

Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)

Director General of Engineering Affairs, MOH

Director General of Royal Hospital

Director General of Khoula Hospital

Director General of Medical Supplies (MOH)

Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)

Hospital Director (Al Nahda Hospital)

Hospital Director (Al Massara Hospital)

The Head of Medical Services in SQU Hospital

The Head of Medical Services in Royal Oman Police

The Head of Medical Services in Ministry of Defence

The Head of Medical Services in The Diwan

The Head of Medical Services in The Sultan's Special Force

The Head of Medical Services in Internal Security Services

The Head of Medical Services in Petroleum Development of Oman

The Head of Medical Services in LNG Oman

**ALL PRIVATE PHARMACIES & DRUG STORES**

After Compliments,

Please find attached our Circular No. **55**..... dated **03/03/20** Regarding Safety Alerts of Anesthesia Systems of (Mfr: GE Healthcare)

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information

# Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs  
and Drug Control

MUSCAT



سلطنة عمان  
وزارة الصحة  
المديرية العامة للصيدلانية  
والرقابة الدوائية  
مسقط

Circular No. 55 / 2020

08 -07-1441 H

03 -03-2020

Ref: 36/2020

## Safety Alerts of Anesthesia Systems (Carestation 620/650/650c (A1 & A2)) from GE Healthcare

Source of Recall	Gulf Health Council																								
Product	Anesthesia Systems (Carestation 620/650/650c (A1 & A2))																								
Manufacturer	GE Healthcare																								
Local Agent	Muscat Pharmacy & Store L.L.C																								
The affected products	<table border="1"><thead><tr><th colspan="3">Affected Devices - WU Manufactured</th></tr><tr><th>Year (YY)</th><th>Fiscal Week (FW)</th><th>Manufacture Site (SA)</th></tr></thead><tbody><tr><td>2018</td><td>34 to 52</td><td>WA</td></tr><tr><td>2019</td><td>01 to 24</td><td>WA</td></tr><tr><th colspan="3">Affected Devices - MA Manufactured</th></tr><tr><th>Year (YY)</th><th>Fiscal Week (FW)</th><th>Manufacture Site (SA)</th></tr><tr><td>2018</td><td>34 to 52</td><td>MA</td></tr><tr><td>2019</td><td>01 to 30</td><td>MA</td></tr></tbody></table>	Affected Devices - WU Manufactured			Year (YY)	Fiscal Week (FW)	Manufacture Site (SA)	2018	34 to 52	WA	2019	01 to 24	WA	Affected Devices - MA Manufactured			Year (YY)	Fiscal Week (FW)	Manufacture Site (SA)	2018	34 to 52	MA	2019	01 to 30	MA
Affected Devices - WU Manufactured																									
Year (YY)	Fiscal Week (FW)	Manufacture Site (SA)																							
2018	34 to 52	WA																							
2019	01 to 24	WA																							
Affected Devices - MA Manufactured																									
Year (YY)	Fiscal Week (FW)	Manufacture Site (SA)																							
2018	34 to 52	MA																							
2019	01 to 30	MA																							
Reason for Recall	The Manufacturer has become aware that there is a potential for a loose cable connection inside specific manufactured anaesthesia devices. This would cause a loss of mechanical ventilation and the system will provide high priority audio and visual alarms. Loss of mechanical ventilation could lead to hypoxia if the clinician does not intervene																								
Action	<ol style="list-style-type: none"><li>1. Kindly check your stock, contact your local agent for remedial action</li><li>2. . If you observe the message – “Ventilate Manually!” change from mechanical to manual ventilation. At any time, the clinician may use a self-inflating bag to ventilate the patient and/or switch to another anaesthesia device.</li><li>3. Perform the planned maintenance (PM) every 12- months at a minimum per the User’s Reference Manual, which includes inspection of the cable connection.</li></ol>																								
Product image																									
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical Device to Director of Medical Device Control contact E-mail <a href="mailto:dg-padc@moh.gov.om">dg-padc@moh.gov.om</a>																								

Directorate General of Pharmaceutical Affairs & Drug Control  
Sultanate of Oman

Dr. Mohammed Hamdan Al Rubaie

DIRECTOR GENERAL

