


GP Clinic



Doctor's room 3.0m X 3.0m	Desk Examination table Medical curtain Wash basin
Reception/ waiting area	Desk Appropriate waiting space for patients
Treatment room 3.0m X 3.0m	Treatment couch Medical curtain Crash cart Wash basin
Observation room 4.0m X 3.0m	Accommodate 2 beds separated by curtains Nurse station Wash basin
Store	For medical consumables
Dirty utility	Wash basin Sink for clinging equipment
Sterilization room	Two sinks for soaking / rinsing equipment Autoclave Class B Wash basin
Medical waste room	With secure lock for storing medical waste
X ray room	Please refer to MOH web site for requirements and specifications
laboratory	25square meter for Hematology and biochemistry machines only, If additional equipment is needed, a larger area must be provided (<i>Please refer to MOH web site</i>)

